



Recombinant *Saccharomyces cerevisiae* Transposon Ty1-DR3 Gag-Pol polyprotein (TY1B-DR3), partial

Product Code	CSB-MP859919SVG
Abbreviation	TY1B-DR3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99231
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Source	Mammalian cell
Target Names	TY1B-DR3
Protein Names	Recommended name: Transposon Ty1-DR3 Gag-Pol polyprotein Alternative name(s): Gag-Pol-p199 TY1A-TY1B Transposon Ty1 TYA-TYB polyprotein p190 Cleaved into the following 4 chains: 1. Capsid protein Short name= 2. CA Alternativ
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.