



Recombinant Human Homeobox protein Meis3 (MEIS3)

Product Code	CSB-MP859938HU
Abbreviation	MEIS3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99687
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MARRYDELPH YPGIVDGPAA LASFPETVPA VPGPYGPHRP PQPLPPGLDS DGLKREKDEI YGHPLFPLLA LVFEKCELAT CSPRDGAGAG LGTPPGGDVC SSDSFNEDIA AFAKQVRSER PLFSSNPELD NLMIQAIQVL RFHLLLELEKV HDLCDNFCHR YITCLKGKMP IDLVIEDRDG GCREFEDYYP ASCPSLPDQN NMWIRDHEDS GSVHLGTPGP SSGGLASQSG DNSSDQGDGL DTSVASPSSG GEDEDLDQER RRNKKRGIFP KVATNIMRAW LfqHLSHPYP SEEQKKQLAQ DTGLTILQVN NWFINARRRI VQPMIDQSNR TGQGAAFSPE GQPIGGYTET QPHVAVRPPG SVGMSLNLEG EWHYL
Source	Mammalian cell
Target Names	MEIS3
Protein Names	Recommended name: Homeobox protein Meis3 Alternative name(s): Meis1-related protein 2
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.