



# Recombinant Human Homeobox protein Meis3 (MEIS3)

<b>Product Code</b>	CSB-YP859938HU
<b>Abbreviation</b>	MEIS3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99687
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MARRYDELPH YPGIVDGPAA LASFPETVPA VPGPYGPHRP PQPLPPGLDS DGLKREKDEI YGHPLFPLLA LVFEKCELAT CSPRDGAGAG LGTPPGGDVC SSDSFNEDIA AFAKQVRSER PLFSSNPELD NLMIQAIQVL RFHLLLELEKV HDLCDNFCHR YITCLKGKMP IDLVIEDRDG GCREFEDYYP ASCPSLPDQN NMWIRDHEDS GSVHLGTPGP SSGGLASQSG DNSSDQGDGL DTSVASPSSG GEDEDLDQER RRNKKRGIFP KVATNIMRAW LfqHLSHPYP SEEQKKQLAQ DTGLTILQVN NWFINARRRI VQPMIDQSNR TGQGAAFSPE GQPIGGYTET QPHVAVRPPG SVGMSLNLEG EWHYL
<b>Source</b>	Yeast
<b>Target Names</b>	MEIS3
<b>Protein Names</b>	Recommended name: Homeobox protein Meis3 Alternative name(s): Meis1-related protein 2
<b>Expression Region</b>	1-375
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.