



Recombinant Human Neuronal PAS domain-containing protein 1 (NPAS1)

Product Code	CSB-BP859941HU
Abbreviation	NPAS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99742
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAPYPGSSGG GSEVKCVGGR GASVPWDFLP GLMVKAPSGP CLQAQRKEKS RNAARSRRGK ENLEFFELAK LLPLPGAISS QLDKASIVRL SVTYLRLRRF AALGAPPWGL RAAGPPAGLA PGRRGPAALV SEVFEQHLGG HILQSLDGFV FALNQEGKFL YISETVSIYL GLSQVEMTGS SVFDYIHPGD HSEVLEQLGL RTPTPGPPTP PSVSSSSSSSS SSLADTPEIE ASLTKVPPSS LVQERSFFVR MKSTLTKRGL HVKASGYKVI HVTGRLRAHA LGLVALGHTL PPAPLAELPL HGHMIVFRLS LGLTILACES RVSDHMDLGP SELVGRSCYQ FVHGQDATRI RQSHVDLLDK GQVMTGYRWW LQRAGGFVWL QSVATVAGSG KSPGEHHVLW VSHVLSQAEG GQTPLDAFQL PASVACEEAS SPGPEPEPE PPTGKQAAP AENEAPQTQG KRIKVEPGPR ETKGSEDSGD EDPSSHPATP RPEFTSVIRA GVLKQDPVRP WGLAPPDPP PTLHAGFLP PVVRGLCTPG TIRYGPAELG LVYPHLQRLG PGPALPEAFY PPLGLPYPGP AGTRLPRKGD
Source	Baculovirus
Target Names	NPAS1
Protein Names	Recommended name: Neuronal PAS domain-containing protein 1 Short name= Neuronal PAS1 Alternative name(s): Basic-helix-loop-helix-PAS protein MOP5 Class E basic helix-loop-helix protein 11 Short name= bHLHe11 Member of PAS pro
Expression Region	1-590
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the basic helix-loop-helix (bHLH)-PAS family of transcription factors. Studies of a related mouse gene suggest that it functions in neurons. The exact function of this gene is unclear, but it may play protective or



modulatory roles during late embryogenesis and postnatal development.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.