



Recombinant Mouse GH3 domain-containing protein (Ghdc)

Product Code	CSB-EP859954MO-B
Abbreviation	Ghdc
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99J23
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	IL WQQRSRGARP CWLISLQHRV AWGMLGWAAA WQQWRLDRST LNVGQSQQQA LMWCLKKAQG SCCLPREDTD MRTFRNHLPL TQTSHTQEQE SEETLPSPAS PQYHGDAASLQ ATLLGLITLN KAYPEALAPG STACVTPTSP WPCSVPWLGH ALGRVSPDGA KDPRTLLEA LISPGLRVLE ARTAVELLDV FVGLEADGEE LAEAIAGIL GTLLPKRAAE LKEALEQGPR GLARRLWPKL QVVVTLDSGG QAEVAALRV LWCQGLAFFS PAYAASGGVV ALNLWPERPQ GSYLLPPGVP FIELLPIKEG TQEEAASTLL LTDAQREKEY ELVLNHTSL TRCRLGDVVQ VVGTYNQCPV VRFTCRLGQT LNVRGEVTDE TVFSVALAQA VGQWPGAKLL DHVCVESRVL DSCEGSAPHY EVFVELRGLR NLSEENRDKL DNCLQEASAQ YKSLRFRGSV GPAKVHLVRP GSFRVLREAL AAFSSSSCRP PEMPRVIRLR HLAQLLQKRV IS
Source	E.coli
Target Names	Ghdc
Protein Names	Recommended name: GH3 domain-containing protein
Expression Region	19-532
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.