



Recombinant Cell division protein FtsZ (ftsZ)

Product Code	CSB-MP860156WAN
Abbreviation	ftsZ
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9ALA3
Product Type	Recombinant Protein
Immunogen Species	Wigglesworthia glossinidia brevipalpis
Purity	>85% (SDS-PAGE)
Sequence	MFPEIELTSD AVIRVVGIGG GGGNAVEHVM RECIEGVDF AVNTDAQALR KTEVSQTVQI GSSITKGLGA GANPEVGKNS AEEDKDALRI ILDGADMVFI ASGMGGGTGT GAAPVIAEIA KDLGILTVAV VTKPFNFEGK KRLIFAEQGI DELSKHVDSL IIPNDKLLK VLGKGISLLD AFSAANDVLK NAVQGI AELI TRPGLINVDF ADVKTMSEM GYAMMGSGIS KGDNR AEESS EIAISSPLLE DIDLSGARGV LVNITAGFDL RLDEF EAVGN KVRFS S DNA TVVIGTSLDP SMNDEL RVTV VATGIGMDKR PDIKLVNST SNKNIMDRFG YRYS DKENTM SKNNNEFSKI KNKTKEDLQS DYLEIPAFLR NQAD
Source	Mammalian cell
Target Names	ftsZ
Protein Names	Recommended name: Cell division protein FtsZ
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>