



Recombinant *Oryza sativa* subsp. *japonica* Putative AP2/ERF and B3 domain-containing protein Os01g0140700 (Os01g0140700, LOC_Os01g04750)

Product Code	CSB-BP860173OFG
Abbreviation	Os01g0140700, LOC_Os01g04750
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9AWS7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEQEAAMVVF SCNSGSGGSS STTDSKQEEE EEEELAAMEE DELIHVVQAA ELRLPSSTTA TRPSSRYKGV VPQPNGRWGA QIYERHARVW LGTFPDEEAA ARAYDVAALR FRGRDAVTNR APAAEGASAG ELAFLAAHSK AEVVDMLRKH TYDDELQQGL RRGSRAPTP RWAREPLFEK AVTPSDVGKL NRLVVPKQQA ERHFPFLRR HSSDAAGKGV LLNFEDGDGK VWRFRYSYWN SSQSYVLTKG WSRFVREKGL RRGDTVAFSR SAAAWGTEKH LLIDCKKMER NNLATVDDDA RVVVKLFGVD IAGDKTR
Source	Baculovirus
Target Names	Os01g0140700
Protein Names	Recommended name: Putative AP2/ERF and B3 domain-containing protein Os01g0140700
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.