



Recombinant Enolase 2 (ENO2)

| | |
|--------------------------|--|
| Product Code | CSB-MP860239TOV |
| Abbreviation | ENO2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9BPL7 |
| Product Type | Recombinant Protein |
| Immunogen Species | Toxoplasma gondii |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MVAIKDITAR QILDSRGNPT VEVDLLTDGG CFRAAVPSGA STGIYEALEL RDKDQTKFMG KGVMKA VENI HKIIKPALIG KDPCDQKGID KLMVEELDGT KNEWGWCKSK LGANAILAVS MACCRAGAAA KGMPLYKYIA TLAGNPTDKM VMPVPPFNVI NGGSHAGNKV AMQEFMIAPV GASTIQEAIQ IGAEVYQHLK VVIKKKYGLD ATNVGDEGGF APNISGATEA LDLLMEAIKV SGHEGKVKIA ADVAASEFFL QDDKVYDLDF KTPNNDKSQR KTGEELRNLY KDLCQKYPFV SIEDPFDQDD FHSYAQLTNE VGEKVQIVGD DLLVTNPTRI EKAVQEKACN GLLLKVNQIG TVSESIACQ LAQKNKWGVM VSHRSGETED SFIADLVVGL RTGQIKTGAP CRSERLCKYN QLMRIEESLG SDCQYAGAGF RHPN |
| Source | Mammalian cell |
| Target Names | ENO2 |
| Protein Names | Recommended name: Enolase 2 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 2 2-phosphoglycerate dehydratase 2 |
| Expression Region | 1-444 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |