



# Recombinant Caiman crocodilus apaporiensis L-lactate dehydrogenase B chain (LDHB)

<b>Product Code</b>	CSB-MP860313CMAK
<b>Abbreviation</b>	LDHB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q98SL1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caiman crocodilus apaporiensis (Rio Apaporis caiman)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ATLKEKLIT HIASESTIPN NKVTVVGVGQ VGMACAVSIL GKSLCDELAL VDVLEDKLGK EMMDLQHGS LFLQTHKIVAD KDYAVTANSK IVVVTAGVRQ QEGESRLNLV QRNVNVFKFI IPQIIKYSPD CILVVSNPV DILTYVTWKL SGLPKHRVIG SGCNLD SARF RYLMSEK LGI HPSSCHGWIL GEHGDSSVAV WSGVNVAGVS LQELNPAMGT DRDSEKWKEV HKQVVESAYE VIKLGKGTNW AIGFSVADLL ETMMKNLNRV QPVSTMVKGM YGIENEVFLS LPCVLSASGL TSVINQKLKD DEVAQLRSSA DTLWSIQKDL KDL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	LDHB
<b>Protein Names</b>	Recommended name: L-lactate dehydrogenase B chain Short name= LDH-B EC= 1.1.1.27
<b>Expression Region</b>	2-333
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.