



Recombinant *Saccharomyces cerevisiae* 3,4-dihydroxy-2-butanone 4-phosphate synthase (RIB3)

Product Code	CSB-YP860330SVG
Abbreviation	RIB3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99258
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MFTPIDQAIE HFKQNKVIV MDDAGRENEG DLICAAENV S TEQMAFLVRH SSGYVCAPMT NAIADKLDLP LLRTGMKFES NDDDRHGTAY TITVDVAQGT TTGISAHDRS MTCRALADSS STPKSFLKPG HICPLRAADG GVLQRRGHT E AGVDLCKLSG LSPVAVIGEL VNDDEQGTMM RLNDCCQAFGK KHGIPLISIE ELAQYLKK
Source	Yeast
Target Names	RIB3
Protein Names	Recommended name: 3,4-dihydroxy-2-butanone 4-phosphate synthase Short name= DHBP synthase EC= 4.1.99.12
Expression Region	1-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.