



# Recombinant Human TSC22 domain family protein 3 (TSC22D3)

<b>Product Code</b>	CSB-BP860343HU
<b>Abbreviation</b>	TSC22D3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99576
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNTEMYQTPM EVAVYQLHNF SISFFSSLLG GDVVSVKLDN SASGASVVAI DNKIEQAMD L VKNHLMYAVR EEVEILKEQI RELVEKNSQL ERENTLLKTL ASPEQLEKFQ SCLSPEEPAP ESPQVPEAPG GSAV
<b>Source</b>	Baculovirus
<b>Target Names</b>	TSC22D3
<b>Protein Names</b>	Recommended name: TSC22 domain family protein 3 Alternative name(s): DSIP-immunoreactive peptide Short name= Protein DIP Short name= hDIP Delta sleep-inducing peptide immunoreactor Glucocorticoid-induced leucine zipper protei
<b>Expression Region</b>	1-134
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein shares significant sequence identity with the murine TSC-22 and Drosophila shs, both of which are leucine zipper proteins, that function as transcriptional regulators. The expression of this gene is stimulated by glucocorticoids and interleukin 10, and it appears to play a key role in the anti-inflammatory and immunosuppressive effects of this steroid and chemokine. Transcript variants encoding different isoforms have been identified for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.