



# Recombinant Human Pituitary homeobox 2 (PITX2)

<b>Product Code</b>	CSB-YP860349HU
<b>Abbreviation</b>	PITX2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99697
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	METNCRKLVS ACVQLGVQPA AVECLFSKDS EIKKVEFTDS PESRKEAASS KFFPRQHPGA NEKDKSQQ GK NEDVGAEDPS KKKRQRRQRT HFTSQQLQEL EATFQRNRY P DMSTREEIAV WTNLTEARVR VWFKNRRAKW RKRERNQQA E LCKNGFGPQF NGLMQPYDDM YPGYSYNNWA AKGLTSASLS TKSFPFFNSM NVNPLSSQSM FSPPNSSISM SMSSSMVPSA VTGVPGSSLN SLNNLNNLSS PSLNSAVPTP ACPYAPPTPP YVYRDTCNSS LASLRLKAKQ HSSFGYASVQ NPASNLSACQ YAVDRPV
<b>Source</b>	Yeast
<b>Target Names</b>	PITX2
<b>Protein Names</b>	Recommended name: Pituitary homeobox 2 Alternative name(s): ALL1-responsive protein ARP1 Homeobox protein PITX2 Paired-like homeodomain transcription factor 2 RIEG bicoid-related homeobox transcription factor Solurshin
<b>Expression Region</b>	1-317
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the RIEG/PITX homeobox family, which is in the bicoid class of homeodomain proteins. The encoded protein acts as a transcription factor and regulates procollagen lysyl hydroxylase gene expression. This protein plays a role in the terminal differentiation of somatotroph and lactotroph cell phenotypes, is involved in the development of the eye, tooth and abdominal organs, and acts as a transcriptional regulator involved in basal and hormone-regulated activity of prolactin. Mutations in this gene are associated with Axenfeld-Rieger syndrome, iridogoniodysgenesis syndrome, and sporadic cases of Peters anomaly. A similar protein in other



vertebrates is involved in the determination of left-right asymmetry during development. Alternatively spliced transcript variants encoding distinct isoforms have been described.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.