



# Recombinant Human Achaete-scute homolog 2 (ASCL2)

<b>Product Code</b>	CSB-MP860358HU
<b>Abbreviation</b>	ASCL2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99929
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDGGTLPRSA PPAPPVPVGC AARRRPASPE LLRCSRRRRP ATAETGGGAA AVARRNERER NRVKLVNLGF QALRQHVPHG GASKKLSKVE TLRSAYEYIR ALQRLLAEHD AVRNALAGGL RPQAVRPSAP RGPPGTTPVA ASPSRASSSP GRGGSSEPGS PRSAYSSDDS GCEGALSPAE RELLDFFSSWL GGY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ASCL2
<b>Protein Names</b>	Recommended name: Achaete-scute homolog 2 Short name= ASH-2 Short name= hASH2 Alternative name(s): Class A basic helix-loop-helix protein 45 Short name= bHLHa45 Mash2
<b>Expression Region</b>	1-193
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene is a member of the basic helix-loop-helix (BHLH) family of transcription factors. It activates transcription by binding to the E box (5 - CANNTG-3 ). Dimerization with other BHLH proteins is required for efficient DNA binding. Involved in the determination of the neuronal precursors in the peripheral nervous system and the central nervous system.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.