



Recombinant Mouse NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 10, mitochondrial (Ndufa10)

Product Code	CSB-MP860381MO
Abbreviation	Ndufa10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99LC3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	LRYGL LAAILGDKTT KKLHEYSRVI TVDGNICSGK NKLAKEIAQQ LGMKHYPEAG IQYSSTTTGD GRPLDIEFSG SCSLEKFYDD PKSNDGNSYR LQSWLYASRL LQYADALEHL LSTGQGVVLE RSIYSDFVFL EAMYNQGYIR KQCVDHYNEI KRLTLPEYLP PHAVIYIDVP VPEVQSRIQK KGDPHEMKVT SAYLQDIENA YKKTFLPKMS EMCEVLVYDS WEAEDPTKVV EDIEYLKYNK GPWLKQDDWT FHYLRMLVQD KTEVLNYTTI PVYLPEITIG AHQGSRIYNS FRELPGRKYA PGYNAEVGDK WIWLK
Source	Mammalian cell
Target Names	Ndufa10
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 10, mitochondrial Alternative name(s): Complex I-42kD Short name= CI-42kD NADH-ubiquinone oxidoreductase 42 kDa subunit
Expression Region	36-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein belongs to the complex I 42kDA subunit family. Mammalian complex I is the first enzyme complex in the electron transport chain of mitochondria. It is composed of 45 different subunits. This protein is a component of the hydrophobic protein fraction and has NADH dehydrogenase activity and oxidoreductase activity. It transfers electrons from NADH to the respiratory chain.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.