




Recombinant Human Catenin beta-1 (CTNNB1)

Product Code	CSB-YP006169HU
Relevance	Key downstream component of the canonical Wnt signaling pathway. In the absence of Wnt, forms a complex with AXIN1, AXIN2, APC, CSNK1A1 and GSK3B that promotes phosphorylation on N-terminal Ser and Thr residues and ubiquitination of CTNNB1 via BTRC and its subsequent degradation by the proteasome. In the presence of Wnt ligand, CTNNB1 is not ubiquitinated and accumulates in the nucleus, where it acts as a coactivator for transcription factors of the TCF/LEF family, leading to activate Wnt responsive genes. Involved in the regulation of cell adhesion. Acts as a negative regulator of centrosome cohesion. Involved in the CDK2/PTPN6/CTNNB1/CEACAM1 pathway of insulin internalization. Blocks anoikis of malignant kidney and intestinal epithelial cells and promotes their anchorage-independent growth by down-regulating DAPK2. Disrupts PML function and PML-NB formation by inhibiting RANBP2-mediated sumoylation of PML . Promotes neurogenesis by maintaining sympathetic neuroblasts within the cell cycle .
Abbreviation	Recombinant Human CTNNB1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35222
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	ATQADLMELDMAMEPDRKAAVSHWQQQSYLDSGIHSGATTTAPSLSGKGNP EEEDVDTSQVLYEWEQGFSSQFTQEQVADIDGQYAMTRAQRVRAAMFPETL DEGMQIPSTQFDAHPTNVQRLAEPSSQMLKHAVVNLINYQDDAELATRAIPEL TKLLNDEDQVVVNKAAMVMVHQLSKKEASRHAIMRSPQMVSIVRTMQNTNDV ETARCTAGTLHNLSHHREGLLAIFKSGGIPALVKMLGSPVDSVLFYAITTLHNL LHQEGAKMAVRLAGGLQKMOVALLNKTNVKFLAITTDCLQILAYGNQESKLILAS GGPQALVNIMRTYTYEKLLWTTSRVLKVLSSVCSSNKPAIVEAGGMQALGLHLT DPSQRLVQNCLWTLRNLSDAATKQEGMEGLLGLTLVQLLGSDDINVVTCAAGIL SNLTCNNYKNKMMVCQVGGIEALVRTVLRAGDREDITEPAICALRHLTSRHQE AEMAQNAVRLHYGLPVVVKLLHPPSHWPLIKATVGLIRNLALCPANHAPLREQ GAIPRLVQLLVRAHQDTQRRTSMGGTQQQFVEGVRMEEIVEGCTGALHILAR DVHNRIVIRGLNTIPLFVQLLYSPIENIQRVAAGVLCELAQDKEAAEAIEAEGATA PLTELLHSRNEGVATYAAAVLFRMSKDPQDYKKRLSVELTSSLFRTEPMAWN ETADLGLDIGAQGEPLGYRQDDPSYRSFHSGGYGQDALGMDPMMHEHEMGG HHPGADYPVDGLPDLGHAQDLMDGLPPGDSNQLAWFDTDL
Research Area	Cancer
Source	Yeast



Target Names	CTNNB1
Protein Names	Beta-catenin
Expression Region	2-781aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	87.4 kDa
Protein Length	Full Length of Mature Protein

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.