



Recombinant Arabidopsis thaliana Peter Pan-like protein (PPAN)

Product Code	CSB-YP860575DOA
Abbreviation	PPAN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ASU7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGRFKNNKKS RVIGKPIAKK NQEDVSHVTG DKIPKSFVFS RMKLAGPVKQ LQMDLRKLML PYTALSLKEK KRNTLRDFLN VSGPMGVTHF LMLSKTASSL SLRVARTPQG PTLTFKIHQY SLASDIAQSQ LRPRCPQDLF KSPPLIVLSG FGSQELHLKL ATIMFQNIFF AIDINTVKLS TCQRLVLLNY NKDTKLIDFR HYSIRLQPVG VSRRIKRVQ NHQVPDLRNL QDVSDVTKA GYGSESEGDE EAATVTLSSD LGRVNKGATK SAVKLQEIGP RMTMQLVKVE EGLCTGGIIF SEYDVDGKKE KLKKKQDEEE EEDSEEEGEE GSEEDGEDMD EDFED
Source	Yeast
Target Names	PPAN
Protein Names	Recommended name: Peter Pan-like protein
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.