



Recombinant *Oryza sativa* subsp. japonica Uroporphyrinogen decarboxylase 1, chloroplastic (Os01g0622300, LOC_Os01g43390)

Product Code	CSB-MP860585OFG
Abbreviation	Os01g0622300, LOC_Os01g43390
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9AXB0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	AARGEDGLPR PPAWMMRQAG RYMAEYQALA KRHPSFRERS ETTTELIVEIT LQPWRAFAPD GVILFSDILT PLPAIGVPFD ISDSKGPVIQ SPVRSEEQVR ELTPIDFEKL RFVGESLKIL RTEIDGQAAL LGFVGAPWTI ATYVVEGGMT NTYTNIKSMC HTAPNVLRLG LSHLADAISE YIIYQVNSGA QCIQIFDSWG GQLPPHVWEQ WSKPYIKQIV NKIKIECPNV PLVLYINGNG GLLERMTDTG VDVIGLDWTV DMADGRRRLG NKISVQGNVD PAFLFSPLPV LTDEIHRVVK SAGPKGHLN LGHGVLVKTP EEAVAHFFDV TRSLRYDTLF QGCVTEVLEP VA
Source	Mammalian cell
Target Names	Os01g0622300
Protein Names	Recommended name: Uroporphyrinogen decarboxylase 1, chloroplastic Short name= UPD1 Short name= URO-D1 EC= 4.1.1.37
Expression Region	51-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.