



# Recombinant Human DnaJ homolog subfamily C member 7 (DNAJC7)

<b>Product Code</b>	CSB-EP860771HU-B
<b>Abbreviation</b>	DNAJC7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99615
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           AAAAECDVV MAATEPELLD DQEAKREAET FKEQGNAYYA KKDYNEAYNY            YTKAIDMCPK NASYYGNRAA TLMMLGRFRE ALGDAQQSVR LDDSFVRGHL            REGKCHLSLG NAMAACRSFQ RALELDHKNA QAQQEFKNAN AVMEYEKIAE            TDFEKRDFRK VVFCMDRALE FAPACHRFKI LKAECLAMLG RYPEAQSVAS            DILRMDSTNA DALYVRGLCL YYEDCIEKAV QFFVQALRMA PDHEKACIAC            RNAKALKAKK EDGNKAFKEG NYKLAYELYT EALGIDPNNI KTNAKLYCNR            GTVNSKLRKL DDAIEDCTNA VKLDDTYIKA YLRRACQYMD TEQYEEAVRD            YEKVYQTEKT KEHKQLLKNA QLELKKSKRK DYYKILGVDK NASEDEIKKA            YRKRALMHHP DRHSGASAEV QKEEEKKFKE VGEAFTILSD PKKKTRYDSG            QDLDEEGMNM GDFDPNNIFK AFFGGPGGFS FEASGPGNFF FQFG         </p>
<b>Source</b>	E.coli
<b>Target Names</b>	DNAJC7
<b>Protein Names</b>	Recommended name: DnaJ homolog subfamily C member 7 Alternative name(s): Tetratricopeptide repeat protein 2 Short name= TPR repeat protein 2
<b>Expression Region</b>	2-494
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.