



Recombinant Human DnaJ homolog subfamily C member 7 (DNAJC7)

Product Code	CSB-MP860771HU
Abbreviation	DNAJC7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99615
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AAAAECDVV MAATEPELLD DQEAKREAET FKEQGNAYYA KKDYNEAYNY YTKAIDMCPK NASYYGNRAA TLMMLGRFRE ALGDAQQSVR LDDSFVRGHL REGKCHLSLG NAMAACRSFQ RALELDHKNA QAQQEFKNAN AVMEYEKIAE TDFEKRDFRK VVFCMDRALE FAPACHRFKI LKAECLAMLG RYPEAQSVAS DILRMDSTNA DALYVRGLCL YYEDCIEKAV QFFVQALRMA PDHEKACIAC RNAKALKAKK EDGNKAFKEG NYKLAYELYT EALGIDPNNI KTNAKLYCNR GTVNSKLRKL DDAIEDCTNA VKLDDTYIKA YLRRACQCYMD TEQYEEAVRD YEKVYQTEKT KEHKQLLKNA QLELKKSKRK DYYKILGVDK NASEDEIKKA YRKRALMHHP DRHSGASAEV QKEEEKKFKE VGEAFTILSD PKKKTRYDSG QDLDEEGMNM GDFDPNNIFK AFFGGPGGFS FEASGPGNFF FQFG
Source	Mammalian cell
Target Names	DNAJC7
Protein Names	Recommended name: DnaJ homolog subfamily C member 7 Alternative name(s): Tetratricopeptide repeat protein 2 Short name= TPR repeat protein 2
Expression Region	2-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.