



Recombinant Human Epidermal growth factor-like protein 8 (EGFL8)

Product Code	CSB-MP860784HU
Abbreviation	EGFL8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99944
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AKGGS LRESQGVCSK QTLVVPLHYN ESYSQPVYKP YLTL CAGRRI CSTYRTMYRV MWREVRREVQ QTHAVCCQGW KKRHPGALTC EAICAKPCLN GGVCVRPDQC ECAPGWGGKH CHVDVDECRT SITLCSHHCF NTAGSFTCGC PHDLVLGVDG RTCMEGSPEP PTSASILSVA VREA EKDERA LKQEIHEL RG RLERLEQWAG QAGAWVRAVL PVPPEELQPE QVAELWGRGD RIESLS DQVL LLEERLGACS CEDNSLGLGV NHR
Source	Mammalian cell
Target Names	EGFL8
Protein Names	Recommended name: Epidermal growth factor-like protein 8 Short name= EGF-like protein 8 Alternative name(s): Vascular endothelial statin-2 Short name= VE-statin-2
Expression Region	26-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.