



Recombinant Mouse Methionine adenosyltransferase 2 subunit beta (Mat2b)

Product Code	CSB-BP860809MO
Abbreviation	Mat2b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99LB6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	VGREKELSI HFVPGCCQLV EEEVNIPSRR VLITGATGLL GRAVYKEFQQ SNWHTVGC GF RRARPKEQV NLLDSEAVHH LIHDFQPHVI VHCAAERRPD VVESQPDAAS QLNVGASGNL AKEAAAIGAF LIYISSDYVF DGTNPPYTEE DIPSPLNLYG KTKLDGEKAV LENNLGAAVL RIPVLYGEVE KLEESAVTVM FDKVQFSNKS ANMDHWQQR F PTHVKDVASV CRQLAEKRML DPSIKGTFHW SGNEQMTKYE MACAIADAFN LPSSHLPIT DSPVIGAQR P KNAQLDCSKL ETLGIGQRTP FRTGIKESLW PFLIDKRWRQ TVFH
Source	Baculovirus
Target Names	Mat2b
Protein Names	Recommended name: Methionine adenosyltransferase 2 subunit beta Alternative name(s): Methionine adenosyltransferase II beta Short name= MAT II beta
Expression Region	2-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein belongs to the methionine adenosyltransferase (MAT) family. MAT catalyzes the biosynthesis of S-adenosylmethionine from methionine and ATP. This protein is the regulatory beta subunit of MAT. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.