



Recombinant Mouse Ras association domain-containing protein 1 (Rassf1)

Product Code	CSB-MP860823MO
Abbreviation	Rassf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99MK9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SAEPELIEL RELAPSGRIG PGRTRLERAN ALRIAPGTTT NPSQQHVPGR GHRFQPAGPT THTWCDLCGD FIWGVVRKGL QCAHCKFTCH YRCRALVCLD CCGPRDLGWD SALERDTNVD EAVERETPD L SQAETE QKIK DYNGQINSNL FMSLNKDGSY TGFIVQLKL VRPVSVSSK KPPSLQDARR GTGRSTAVKR RTSFYLPKDA IKHLHVLRSR RAREVIEALL RKFVVDVDP KRALFETER HGQVYLRKLS DDEQPLKRL LAGPSEKALS FVLKENDSGE VNWDAFSMPE LHNFLRILQR EEEHRLRQIL QKYSRCRQKI QEALHACPLG
Source	Mammalian cell
Target Names	Rassf1
Protein Names	Recommended name: Ras association domain-containing protein 1 Alternative name(s): Protein 123F2
Expression Region	2-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a protein similar to the RAS effector proteins. Loss or altered expression of this gene has been associated with the pathogenesis of a variety of cancers, which suggests the tumor suppressor function of this gene. The inactivation of this gene was found to be correlated with the hypermethylation of its CpG-island promoter region. The encoded protein was found to interact with DNA repair protein XPA. The protein was also shown to inhibit the accumulation of cyclin D1, and thus induce cell cycle arrest. Seven alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.