



# Recombinant Mouse Ubiquilin-4 (Ubqln4)

<b>Product Code</b>	CSB-EP860834MO-B
<b>Abbreviation</b>	Ubqln4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99NB8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAEPSGAETR PQIRVTVKTP KDKEEIVICD QASVKEFKEE ISRRFKAQQD QLVLIFAGKI LKDGDTLSQH GIKDGLTVHL VIKTPQKAQD PVTAAASPPS TPDSASAPST TPASPAAAPV QPCSSGNTTS DAGSGGGPSP VAAEGPSSAT ASILSGFGGI LGLGSLGLGS ANFMELQQQM QRQLMSNPPEM LSQIMENPLV QDMMSNPDL M RHMIMANPQM QQLMERNPEI SHMLNNPELM RQTMELARNP AMMQEMMRNQ DRALSNLESV PGGYNALRRM YTDIQEPMFT AAREQFGNNP FSSLAGNSDN SSSQPLRTEN REPLPNPWSP SPPTSQAPGS GGEGTGGSGT SQVHPTVSNP FGINAASLGS GMFNPEMQA LLQQISENPQ LMQNVISAPY MRTMMQTLAQ NPDFAAQMMV NVPLFAGNPQ LQEQLRLQLP VFLQQMQNPE SLSILTNPRA MQALLQIQQG LQTLQTEAPG LVPSLGSFGT PRTSVPLAGS NSGSSAEAPT SSPGVPATSP PSAGSNAQQQ LMQQMIQLLS GSGNSQVPMP EVRFQQQLEQ LNSMGFINRE ANLQALIATG GDINAAIERL LGSQLS
<b>Source</b>	E.coli
<b>Target Names</b>	Ubqln4
<b>Protein Names</b>	Recommended name: Ubiquilin-4 Alternative name(s): Ataxin-1 interacting ubiquitin-like protein Short name= A1Up Ataxin-1 ubiquitin-like-interacting protein A1U Connexin43-interacting protein of 75 kDa Short name= CIP75
<b>Expression Region</b>	1-596
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.