



Recombinant Mouse UBX domain-containing protein 6 (Ubxn6)

Product Code	CSB-MP860843MO
Abbreviation	Ubxn6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99PL6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MKKFFQEIKA DIKFKSAGPG QKLTDSAGEK TTKGKSPQLA LRQPRQGPTD EAQMAAAAAL ARLEQKQPRA RGPTSQDSIR NQVRKELQAE ATSSNPGAP GTNSVPEPKE EISPHLAVPG VFFICPLTGV TLRRDQRDAH IKQAILSHFS TDPVAASIMK IHTFNDRDR DRVCLGVDTIAK YLDNIHLHPE EEKYQKIKLQ NKVFQERINC LEGSHEFFEA IGFKKVTLPV PDQEGQEEFY VLGEDARAQP QNLARHKQQL LDAEPVRATL DRQLRVFRPS ALASHFELPS DFFSLTAEV KREQRLRTEA VERLSSLRTK AMREKEEQRE LRKYTYALVR VRLPDGCLLQ GTFYAREKLS ALFRFVREAL QNDWLFPFELR ASGGQKLEEN EALALNECGL VPSALLTFSW DASVLEDIRA AGAEPKSVL RPELLAAIEQ LS
Source	Mammalian cell
Target Names	Ubxn6
Protein Names	Recommended name: UBX domain-containing protein 6 Alternative name(s): UBX domain-containing protein 1
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.