



# Recombinant Sheep Keratin, type I cytoskeletal 25 (KRT25)

<b>Product Code</b>	CSB-BP861045SH
<b>Abbreviation</b>	KRT25
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BGM5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Ovis aries (Sheep)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSLRLPSGSR RAGPRPTTGS LRLSGAGASF GAGNACSM PG IGSSFSCAFG SSSSGGNALG GNPCAGFTMN EGLLSGNEK VTMQNLNDRL ASYLENVRAL EEANADLEQK IKGWYEKFGP GSCRGLDHDY SRYLPIEDL KNQIIASTTS NANA VLQIDN ARLTADV FRL KYENELALHQ SVESDVNGLR RVLDVITLCR TDLEIQYETL SEELTYLKKH HKEEMQVLQC AAGGNVNVEM NAAPGVDLTV LLNNMRAEYE ALAEQNRDA EAWFNEKSAS LQQQITEDVG ATTSARNELT EMKRNLQTLE IELQSL LTK HSLECSLTET EGN YCAQLAQ VQAQIGALEE QLHQVRTETE GQKLEYEQLL DIKVHLEKEI ETYCLLIGGD DGACKSGGYK SKDYAAGNMG NQMKDPIRAI VVKVLEEVD QRSKVL TTRL HSLEEK SQSN
<b>Source</b>	Baculovirus
<b>Target Names</b>	KRT25
<b>Protein Names</b>	Recommended name: Keratin, type I cytoskeletal 25 Alternative name(s): Cytokeratin-25 Short name= CK-25 Keratin-25 Short name= K25 Type I inner root sheath-specific keratin-K25irs1 Type I keratin intermediate filament IRS
<b>Expression Region</b>	1-450
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the type I (acidic) keratin family, which belongs to the superfamily of intermediate filament (IF) proteins. Keratins are heteropolymeric structural proteins which form the intermediate filament. These filaments, along with actin microfilaments and microtubules, compose the cytoskeleton of epithelial cells. The type I keratin genes are clustered in a region of chromosome 17q12-q21.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.