



# Recombinant Human RNA-binding protein 4B (RBM4B)

<b>Product Code</b>	CSB-MP861079HU
<b>Abbreviation</b>	RBM4B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BQ04
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVKLFIGNLP REATEQEIRS LFEQYGKVL E CDIKNYGFV HIEDKTAAED AIRNLHHYKL HGVNINVEAS KNKSKASTKL HVGNISPTCT NQELRAKFEE YGPVIECDIV KDYAFVHMER AEDA VEAIRG LDNTEFQGKR MHVQLSTSRL RTAPGMGDQS GCYRCGKEGH WSKECPVDRT GRVADFTEQY NEQYGAVRTP YTMGYGESMY YNDAYGALDY YKRYRVR SYE AVAAAAAASA YNYAEQTMSH LPQVQSTTVT SHLNSTSVDP YDRHLLPNSG AAATSAAMAA AAATTSSYYG RDRSPLRRAA AMLPTVGEGY GYGPESELSQ ASAATRNSLY DMARYEREQY VDRARYSAF
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RBM4B
<b>Protein Names</b>	Recommended name: RNA-binding protein 4B Alternative name(s): RNA-binding motif protein 30 RNA-binding motif protein 4B RNA-binding protein 30
<b>Expression Region</b>	1-359
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.