



Recombinant Human Neuralized-like protein 2 (NEURL2)

Product Code	CSB-EP861118HU-B
Abbreviation	NEURL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BR09
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAASEPVDS GALWGLERPE PPTTRFHRVH GANIRVDPSG TRATRVESFA HGVCFSREPL APGQVFLVEI EEKELGWCGH LRLGLTALDP ASLAPVPEFS LPDLVNLGHT WVFAITRHHN RVPREGRPEA EAAAPSRPPT LLVEPYLRIE QFRIPRDRLV GRSRPGLYSH LLDQLYELNV LPPTARRSRL GVLFCPRPDG TADMHIING EDMGPSARGL PAAQPLYAVV DVFASTKSVR LVQLEYGLPS LQTLCLRLVIQ RSMVHRLAID GLHLPKELKD FCKYE
Source	E.coli
Target Names	NEURL2
Protein Names	Recommended name: Neuralized-like protein 2
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.