



Recombinant Human Coiled-coil domain-containing protein 77 (CCDC77)

Product Code	CSB-BP861119HU
Abbreviation	CCDC77
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BR77
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNFTPTHTPV CRKRTVVSKR GVAVSGPTKR RGMADSLEST PLPSPEDRLA KLHPSKELLE YYQKKMAECE AENEDLLKKL ELYKEACEGQ HKLECDLQQR EEEIAELQKA LSDMQVCLFQ EREHVLRLYS ENDRLRIREL EDKKKIQNLL ALVGTDAGEV TYFCKEPPHK VTILQKTIQA VGECEQSESS AFKADPKISK RRPSRERKES SEHYQRDIQT LILQVEALQA QLGEQTKLSR EQIEGLIEDR RIHLEEIQVQ HQRNQNKIKE LTKNLHHTQE LLYESTKDFL QLRSENQNKE KSWMLEKDNL MSKIKQYRVQ CKKKEDKIGK VLPVMHESHH AQSEYIKSLK DKLVQEKKLS NMYQEQCISL EEELARIREE EGMRRREIFKD RTNKM GKRLQ IMTKRYEALE RRRILEVEGF KTDIKVLRQK LKDLEQMPLYK ATVNARANQD LALLCEVRDS NRRAHKIQGE LKNLKS KVFGL ENELRLC
Source	Baculovirus
Target Names	CCDC77
Protein Names	Recommended name: Coiled-coil domain-containing protein 77
Expression Region	1-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.