



Recombinant Human Zinc finger protein with KRAB and SCAN domains 3 (ZKSCAN3)

Product Code	CSB-EP861122HU-B
Abbreviation	ZKSCAN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BRR0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MARELSESTA LDAQSTEDQM ELLVIKVEEE EAGFPSSPDL GSEGSRRERFR GFRYPEAAGP REALSRLREL CRQWLQPEMH SKEQILELLV LEQFLTILPG NLQSWVREQH PESGEEVVVL LEYLERQLDE PAPQVSGVDQ GQELLCKMA LLTPAPGSQS SQFQLMKALL KHESVGSQPL QDRVLQVPVL AHGGCCREDK VVASRLTPES QGLLKVEDVA LTLTPEWTQQ DSSQGNLCRD EKQENHGSLV SLGDEKQTKS RDLPPAEELP EKEHGKISCH LREDIAQIPT CAEAGEQEGR LQRKQKNATG GRRHICHECG KSFAQSSGLS KHRRIHTGEK PYECEECGKA FIGSSALVIH QRVHTGEKPY ECEECGKAFS HSSDLIKHQR THTGEKPYEC DDCGKTFSQS CSLEHHRIH TGEKPYQCSM CGKAFRRSSH LLRHQRIHTG DKNVQEPEQG EAWKSRMESQ LENVETPMSY KCNECERSFT QNTGLIEHQK IHTGEKPYQC NACGKGFTRI SYLVQHQRSH VGKNILSQ
Source	E.coli
Target Names	ZKSCAN3
Protein Names	Recommended name: Zinc finger protein with KRAB and SCAN domains 3 Alternative name(s): Zinc finger and SCAN domain-containing protein 13 Zinc finger protein 306 Zinc finger protein 309 Zinc finger protein 47 homolog Short name=
Expression Region	1-538
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.