



Recombinant Human DCN1-like protein 5 (DCUN1D5)

Product Code	CSB-EP861134HU-B
Abbreviation	DCUN1D5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BTE7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPVKKKRKSP GVAAVAEDG GLKKCKISSY CRSQPPARLI SGEEHFSSKK CLAWFYEYAG PDEVVGPEGM EKFCEDIGVE PENIIMLVLA WKLEAESMGF FTKEEWLKG M TSLQCDCTEK LQNKFDLRS QLNDISSFKN IYRYAFDFAR DKDQRSLDID TAKSMLALLL GRTWPLFSVF YQYLEQSKYR VMNKDQWYNV LEFSRTVHAD LSNYDEDGAW PVLLDEFVEW QKVRQTS
Source	E.coli
Target Names	DCUN1D5
Protein Names	Recommended name: DCN1-like protein 5 Alternative name(s): DCUN1 domain-containing protein 5 Defective in cullin neddylation protein 1-like protein 5
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.