



Recombinant Human Zinc finger protein GLIS2 (GLIS2)

Product Code	CSB-MP861185HU
Abbreviation	GLIS2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BZE0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MHSLDEPLDL KLSITKLRAA REKRERTLG V VRPRALHREL GLVDDSP TPG SPGSPPSGFL LNSKFPEKVE GRFSAAPLVD LSLSPPSGLD SPNGSSSLSP ERQGN GDLPP VPSASDFQPL RYLDGVPSSF QFFLPLGSGG ALHLPASSFL TPPKDKCLSP DLPLPKQLVC RWAKCNQLFE LLQDLVDHVN DYHVKPEKDA GYCCHWEGCA RHGRGFNARY KMLIHIRTHT NEKPHRCPTC SKSFSRLENL KIHNRSH TGE KPYVCPYEGC NKRYSNSSDR FKHTRTHYVD KPYYCKMPGC HKRYTDPSSL RKHIKAHG HF VSHEQQELLQ LRPPPKPPLP APDGGPYVSG AQIIIPNPAA LFGGPGLPGL PLPLAPGPLD LSALACGNNG GSGGGGGMGP GLPGPVLPLN LAKNPLLPSP FGAGGLGLPV VSLLAGAAGG KAEGEKGRGS VPTRALGMEG HKTPLERTES SCSRSPDGL PLLPGTVLDL STGVNSAASS PEALAPGWVV IPPGSVLLKP AVVN
Source	Mammalian cell
Target Names	GLIS2
Protein Names	Recommended name: Zinc finger protein GLIS2 Alternative name(s): GLI-similar 2 Neuronal Krueppel-like protein
Expression Region	1-524
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.