



Recombinant *Saccharomyces pastorianus* Adenylyl-sulfate kinase (MET14)

| | |
|--------------------------|---|
| Product Code | CSB-YP861216SAJ |
| Abbreviation | MET14 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9C2Y6 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Saccharomyces pastorianus</i> (Lager yeast) (<i>Saccharomyces cerevisiae</i> x <i>Saccharomyces eubayanus</i>) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MATNITWHPN LTYDERKELR KQDGCTVWLT GLSASGKSTI ACALEQLLLQ KNLSAYRLDG DNIRFGLNKD LGFSEKDRNE NIRRISEVSK LFADSCAVSI TSFISPYRVD RDRARDLHKE AGLKFIEIFV DVPLEVAEQR DPKGLYKKAR EGVIKEFTGI SAPYEAPKAP ELHLRTDQKT VEECAAIYE YLVNEKIIRK HL |
| Source | Yeast |
| Target Names | MET14 |
| Protein Names | Recommended name: Adenylyl-sulfate kinase EC= 2.7.1.25 Alternative name(s): ATP adenosine-5'-phosphosulfate 3'-phosphotransferase Adenosine-5'-phosphosulfate kinase Short name= APS kinase |
| Expression Region | 1-202 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |