



Recombinant Arabidopsis thaliana DNA replication complex GINS protein PSF2 (GINS2)

Product Code	CSB-EP861241DOA
Abbreviation	GINS2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C7A8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAGQTDPHIS LFSPQEVEFM AEDELVEIVP NMNMEQLNFI SGDFGRFIPQ IPTKVPLWLA VALKRRGKCT FRPPGWMSVD NLTQILEAER ESQSTFQALP FSYVEIARLL FDHARDDIPD MYMVRSLVED IRDVRLHKLE TNLGSFQGTS AVKISNVSAM EVNIVRPFVI RALEAFYKHD KPEADVDRDT RSSRQQRETN NEPRRPLRQR
Source	E.coli
Target Names	GINS2
Protein Names	Recommended name: DNA replication complex GINS protein PSF2
Expression Region	1-210
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.