



Recombinant Arabidopsis thaliana Hydroxyacylglutathione hydrolase 3, mitochondrial (GLY3)

Product Code	CSB-EP861252DOA-B
Abbreviation	GLY3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C8L4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	LFRQLFEN ESSTFTYLLA DVSHDPKPAL LIDPVDKTV D RDLKLIDELG LKLIYAMNTH VHADHVTGTG LLKTKLPGVK SVISKASGSK ADLFLEPGDK VSIQDIYLEV RATPGHTAGC VTYVTGEGAD QPQPRMAFTG DAVLIRGCGR TDFQEGSSDQ LYESVHSQIF TLPKDTLIYP AHDYKGFVVS TVGEE MQHNP RLTKDKETFK TIMSNLNLSY PKMIDVAVPA NMVCGLQDVP SQAN
Source	E.coli
Target Names	GLY3
Protein Names	Recommended name: Hydroxyacylglutathione hydrolase 3, mitochondrial EC=3.1.2.6 Alternative name(s): Glyoxalase II Short name= Glx II
Expression Region	53-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.