



Recombinant Arabidopsis thaliana Probable histone chaperone ASF1A (ASF1A)

Product Code	CSB-EP861262DOA-B
Abbreviation	ASF1A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C9M6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSAIKITNVA VLHNPAPFVS PFQFEISYEC LNSLKDDLEW KLIYVGS AED ETYDQLLESV LVGPVNVGNY RFV FQADPPD PSKIQEEDII GVTVLLLTCS YMGQEFLRVG YYVNNDYEDE QLKEEPPTKV LIDKVQRN IL SDKPRVTKFP IDFHPEEEQT AATAAPPEQS DEQQPNVNGE AQVLPDQSVE PKPEES
Source	E.coli
Target Names	ASF1A
Protein Names	Recommended name: Probable histone chaperone ASF1A Alternative name(s): Anti-silencing function protein 1-like protein a Short name= Anti-silencing function 1a protein S-locus protein 7 Short name= AtSP7 Silencing group A pro
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.