



Recombinant Mouse Protein FAM212A (Fam212a)

Product Code	CSB-MP861480MO
Abbreviation	Fam212a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9CX62
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MHSARLDSFL SQLRWELLCA RDTGSPMSG PLQPKPRTDQ NVQPKRQFRA SDVLEEDSVC CVEEEEEEEGL VAEDKGLPLG CPREHALDWD SGFSEVSGST WREEEPSVPQ RQAPRERPPH SQRFSVSDIP MRSRAAVTNI PPAHRPRPKS TPDACLEHWQ GLEAEDWTAA LLNRGRSRQP LVLGDNCFAD LVHNWMEPE ATSEGS DGDV PRARARPPQF LLGLSEQLRR RLARARRTAM ASKRLSCPPR SEPDLPADIS RFAALMNCRS RQPIIYNDVS YL
Source	Mammalian cell
Target Names	Fam212a
Protein Names	Recommended name: Protein FAM212A
Expression Region	1-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.