



Recombinant Mouse Protein lin-37 homolog (Lin37)

Product Code	CSB-BP861574MO
Abbreviation	Lin37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9D8N6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MFPVKVKVEK SEMEMAKARN QLDAVLQCLL EKSHMDRERL DEEAGKTPLD THNKDCSIAA TGKRPSARFP HQRRKKRREM DDGLAEGGPQ RSNTYVIKLF DRSVDLAQFS ENTPLYPICR AWMRNSPTVR ERERSPGSPL PPLPEDGEGS EVINSKNRDV YKLPPPTAPG PLGDACRSRI PSPLQPETEG TPDDEPSEPE PSPSTLIYRN MQRWKRIRQR WKEASHRNQL RYSESMKILR EMYDRQ
Source	Baculovirus
Target Names	Lin37
Protein Names	Recommended name: Protein lin-37 homolog Alternative name(s): Antolefinin Antolefinine
Expression Region	1-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein expressed in the eye.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.