



Recombinant Mouse General transcription factor IIH subunit 1 (Gtf2h1)

Product Code	CSB-EP861603MO
Abbreviation	Gtf2h1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DBA9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MATSSEEVLL IVKKVRQKKQ DGALYLMAER IAWAPEGKDR FTISHMYADI KCQKISPEGK AKIQLQLVLH AGDTTNFHS NESTAVKERD AVKDLLQQLL PKFKRKANKE LEEKNRMLQE DPVLFQLYKD LVVSQVISAE EFWANRLNVN ATDSSTSSHK QDVGISAAFL ADVRPQTDGC NGLRYNLTS D IIESIFRTYP AVKMKYAETV PHNMTEKEFW TRFFQSHYFH RDRLNTGSKD LFAECAKIDE KGLKTMVSLG VKNPMLDLTS LEDKPLDEGY SISSVPSTSN SKSIKENSNA AIIKRFNHHS AMVLAAGLRK QQAQNGQNGE PSSVDGNSGD TDCFQPAVKR AKLQESIEYE DLGNNNSVKT IALNLKKS DR YYHGPTPIQS LQYATSQDII NSFQSIRQEM EAYTPKLTQV LSSSAASSTI TALSPGGALM QGGTQQAVNQ MVPNDIQSEL KHLYVAVGEL LRHFWSFCFPV NTPFLEEKVV KMKSNLRFQ VTKLCPFQEK IRRQYLSTNL VSHIEEMLQT AYNKLHTWQS RRLMKKT
Source	E.coli
Target Names	Gtf2h1
Protein Names	Recommended name: General transcription factor IIH subunit 1 Alternative name(s): Basic transcription factor 2 62 kDa subunit Short name= BTF2 p62 General transcription factor IIH polypeptide 1 TFIIH basal transcription factor compl
Expression Region	1-547
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.