



# Recombinant *Xenopus laevis* Calcium/calmodulin-dependent protein kinase type II delta chain (camk2d)

<b>Product Code</b>	CSB-BP861635XBE
<b>Abbreviation</b>	camk2d
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9DG02
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus laevis</i> (African clawed frog)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASTTCTRFT DEYQLFEELG KGAFSVVRRRC IKINIGQEYA AKIINTKKLS ARDHQKLERE AKICRLLKHP NIVRLHDSIS EEGFHYLVFD LVTGGELFED IVAREYYSEA DASHCIQQIL ESVNHCHLNG IVHRDLKPEN LLLASKLKGA AVKLADFGLA IEVQGDQQAW FGFAGTPGYL SPEVLRKDPY GKPVDMWACG VILYILLVGY PPFWDEDQHR LYQQIKAGAY DFPSPWDTV TPEAKDLINK MLTINPAKRI NATEALRHPW ICQRSTVASM MHRQETVDCL KKFNARRKLLK GAILTTMLAT RNFSAKSLLK KPDGVKESTE SSNTTIEDDED VKARKQEIIK VTEQLIEAIN NGDFEAYTKI CDPGLTSFEP EALGNLVEGM DFHRFYFENA LSKSNKPVHT IILNPHVHLI GEDAACIAYI RLTQYLDSDAG MPKTMQSEET RVWHRRDGKW QNVHFHRSGS PTIPN
<b>Source</b>	Baculovirus
<b>Target Names</b>	camk2d
<b>Protein Names</b>	Recommended name: Calcium/calmodulin-dependent protein kinase type II delta chain Short name= CaM kinase II subunit delta Short name= CaM-kinase II delta chain Short name= CaMK-II subunit delta EC= 2.7.11.17
<b>Expression Region</b>	1-475
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .