



Recombinant Human hepatitis A virus genotype IIIA Genome polyprotein

Product Code	CSB-EP861650HFAY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DWR1
Product Type	Recombinant Protein
Immunogen Species	Human hepatitis A virus genotype IIIA (isolate NOR-21) (HHAV) (Human hepatitis A virus (isolate Human/Norway/NOR-21/1998))
Purity	>85% (SDS-PAGE)
Sequence	MNMSRQGIFQ TVGSGLDHIL SLADVEEEQM IQSVDRTAVT GASYFTSVDQ SSVHTAEVGS HQPEPLKTSV DKPGSKRTQG EKFFLIHSAD WLTTALFHE VAKLDVVKLL YNEQFAVQGL LRYHTYARFG IEIQVQINPT PFQQGGLICA MVPGDQSYGS IASLTVYPHG LLNCNINNVV RIKVPFIYTR GAYHFKDPQY PVWELTIRVW SELNIGTGTS AYTSLNLVLAR FTDLELHGLT PLSTQ
Source	E.coli
Protein Names	Recommended name: Genome polyprotein Cleaved into the following 18 chains: 1. Protein VP0 Alternative name(s): VP4-VP2 Protein VP4 Alternative name(s): P1A Virion protein 4 Protein VP2 Alternative name(s): P1B Virion pr
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.