



Recombinant Mesocricetus auratus Interleukin-17B (IL17B)

Product Code	CSB-BP861679MRG
Abbreviation	IL17B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9EQI6
Product Type	Recombinant Protein
Immunogen Species	Mesocricetus auratus (Golden hamster)
Purity	>85% (SDS-PAGE)
Sequence	RNTKGKRK GQVRPGPLAP GPHQVPLDLV SRVKPYARME EYERNLGEMV AQLRNSSEPA KRRCEVNLQL WLSNKRSLSP WGYSINHDP S RIPADLPEAR CLCLGCVNPF TMQEDRSMVS VPVFSQVPVR RRLCPPPPRP GPCRHRVVME TIAVGCTC
Source	Baculovirus
Target Names	IL17B
Protein Names	Recommended name: Interleukin-17B Short name= IL-17B Alternative name(s): Neuronal interleukin-17-related factor
Expression Region	23-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.