



Recombinant Arabidopsis thaliana Germin-like protein subfamily 1 member 14 (At5g39110)

Product Code	CSB-MP861800DOA
Abbreviation	At5g39110
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FID0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	NDPSPLQD FCVAIGDLKN GVFVNGKFCK DPKQAKAEDF FYSGLNQAGT TNNKVKSNTV TVNVDQIPGL NTLGISLVRI DYAPYGQNP HTHPRATEIL VLVEGTLYVG FVSSNQDNNR LFAKVLNPGD VFVFPIMIH FQVNIKTPA VAFAGLSSQN AGVITIADTV FGSTPPINPD ILAQAFQLDV NNVKDLEAKF KN
Source	Mammalian cell
Target Names	At5g39110
Protein Names	Recommended name: Germin-like protein subfamily 1 member 14
Expression Region	23-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.