



Recombinant Arabidopsis thaliana F-box protein PP2-A14 (PP2A14)

Product Code	CSB-BP861811DOA
Abbreviation	PP2A14
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FJ80
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGAASSSIVR SEPFAGKLCG LEDVPENCIT AMFMYMEPPE ICLLARVNKS FHRASRSDAV WEDKLPSTNYK FLVRRILEDQ QQVGVKDKLI YRKKEIYARL CRPNLFDTGT KEAWLDRSG KVFLAISPKA MKITGIDDRR YWEHISSEDES RFGSITYLRQ IWWLEAVGKI RFEFAPGKYS LLFKIQLGKP IRKCGRKTCS LDQVHGWDIK PVRFQLSTSD GQCAMSERHL DESGRWVYHH AGDFVVENQN SPVWVKFSML QIDCTHTKGG LCLDCVIICP FEYRGKYKYS D
Source	Baculovirus
Target Names	PP2A14
Protein Names	Recommended name: F-box protein PP2-A14 Alternative name(s): Protein PHLOEM PROTEIN 2-LIKE A14 Short name= AtPP2-A14 SKP1-interacting partner 13
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.