



Recombinant Arabidopsis thaliana Mitochondrial outer membrane protein porin 2 (VDAC2)

Product Code	CSB-YP861818DOA
Abbreviation	VDAC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FJX3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	SKGPGGLFTD IGKKAKDLLT RDYNSDQKFS ISTYSASGVA LTSTALKKGG VHAADVATQY KYKNALFDVK IDTDSSVLTT VTLTEILPST KAIASFKVPD YNSAKLEVQY FHDHATVTAA AALKQNPLID ITATLGSPVI SFGAEAGYDT TSKTFTKYNA GISVTKPDAC LSIILGDKGD SLKASYLHHF DEFKRRTAAVG EVYRKfstNE NTITVGGLYA IDHSTAVKAK LNNHGTALGAL LQHEVLPRSL VTVSSEIDTK ALEKHPRFGL SLALKP
Source	Yeast
Target Names	VDAC2
Protein Names	Recommended name: Mitochondrial outer membrane protein porin 2 Alternative name(s): Voltage-dependent anion-selective channel protein 2 Short name= AtVDAC2 Short name= VDAC-2
Expression Region	2-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.