



Recombinant Arabidopsis thaliana ZF-HD homeobox protein At5g65410 (At5g65410)

Product Code	CSB-BP861828DOA
Abbreviation	At5g65410
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FKP8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEFEDNNNNN DEEQEEDMNL HEEEEEDDDAV YDSPPLSRVL PKASTESHET TGTTSTGGGG GFMVVHGGGG SRFRFRECLK NQAVNIGGHA VDGCGEFMPA GIEGTIDALK CAACGCHRNH HRKELPYFHH APPQHQP PPP PPGFYRLPAP VSYP PPSQA PPLQLALPPP QRERSEDPME TSSAEAGGGI RKRHRKFTA EQKERMLALA ERIGWRIQRQ DDEVIQRFCQ ETGVPRQVLK VWLHNNKHTL GKSP SPLHHH QAPPPPPQS SFHHEQDQP
Source	Baculovirus
Target Names	ZHD1
Protein Names	Recommended name: ZF-HD homeobox protein At5g65410
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.