



Recombinant Arabidopsis thaliana Putative F-box/FBD/LRR-repeat protein At5g44960 (At5g44960)

Product Code	CSB-YP861835DOA
Abbreviation	At5g44960
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FLA1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEECDYINEL PDSLLTQILL DLRTKDSVKT SVSSKRWRNL WLNVPGLDLF SLQFTNPHE EGLIKFMDRF MESNCRSRLQ KFMIRYFECN GYRDRFMELI GTVVDGCIQH LYVYMHTCNR VDFIRQNIYK SKTLVSLKLY NVELKNPDFV VSLPCLKILK LMKICYGEDG PLVVEKLISG CPVLEDLELI KPFDILTQDV ILFLRVSSQT LKSLRLYFAT NGSVEIDAPR LKYMTFYESR FDRIMVKNMS SLFSIEIRAK SSFEYGGLLK AEDPRKRAII CDFLTVISSV RHMIISGSIL EELHSYSKLG WIPQFRNLYH LQASFFGTSL QLLPTFLESC PNLKNLIMDY GAFKEENIDF HEVPQCLIST LEYVHINKLM MMEQSGIKLV NYFIENSAVL KKLTLRFSFF SSIESESYKK LLTSTKLSPT CQVIFV
Source	Yeast
Target Names	At5g44960
Protein Names	Recommended name: Putative F-box/FBD/LRR-repeat protein At5g44960
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.