



Recombinant Arabidopsis thaliana Peroxidase 52 (PER52)

Product Code	CSB-BP861836DOA
Abbreviation	PER52
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FLC0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	DNNYVVEAQL TTNFYSTSCP NLLSTVQTAV KSAVNSEARM GASILRLFFH DCFVNGCDGS ILLDDTSSFT GEQNAAPNRN SARGFNVIDN IKSAVEKACP GVVSCADILA IAARDSVVAL GGPWNVVKVG RRDARTASQA AANSNIPAPT SSLSQLISSF SAVGLSTRDM VALSGAHTIG QSRCTNFRAR IYNETNINAA FATTRQRTCP RASGSGDGNL APLDVTTAAS FDNFYFKNLM TQRGLLHSDQ VLFNGGSTDS IVRGYSNNPS SFNSDFTAAM IKMGDISPLT GSSGEIRKVC GRTN
Source	Baculovirus
Target Names	PER52
Protein Names	Recommended name: Peroxidase 52 Short name= Atperox P52 EC= 1.11.1.7 Alternative name(s): ATP49
Expression Region	21-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.