



Recombinant Arabidopsis thaliana F-box/FBD/LRR-repeat protein At5g22700 (At5g22700)

Product Code	CSB-BP861862DOA
Abbreviation	At5g22700
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FNJ1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSNRGDRISS LPDELLCQIL SNLPTKNAVT TSILSTRWRS IWLSTPVLDI DIDAFDDATT FVSFASRFLE FSKDSCLHKF KLSVERDDVD MCTIMPWIQD AVNRRIQHLE VDCRFSHFHE AVYLTLYLSE TLVSLRLHFV TLHRYEFVSL PNLKVMHLEE NIYYCLETLE NFISSCPVLE DLTVVRIVDI ITEKILRVRS RSLNSLKLVL DSSNGWFIDD IDEWKVIIDA PRLAYLSLKD DQSASFVSN LGSSAKADIK VSFNVNDIWD LPVTFERSNV GKLLTGLSSI RDLTISGTTL MIICHYLKHE PMPQFCNMVR LNARFYDCDL EMLPCVLESC PNLKSLVLGL LSSAETEQQS RVSSVPPCFL SSLEFVEIRS RLCRKRYVMK VARYFAKNSV MLKKFVYVGR VSIQEEDLLA LLWRYSICEI EVRGLRR
Source	Baculovirus
Target Names	At5g22700
Protein Names	Recommended name: F-box/FBD/LRR-repeat protein At5g22700
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.