



Recombinant Arabidopsis thaliana Transcription factor MYB113 (MYB113)

Product Code	CSB-EP861865DOA-B
Abbreviation	MYB113
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FNV9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGESPKGLRK GTWTTEEDIL LRQCIDKYGE GKWHRVPLRT GLNRCRKSCR LRWLNLYLKPS IKRGKLCSD EVDLVLRLHKL LGNRWSLIAG RLPGRRTANDV KNYWNTLSK KHDERCCKTK MINKNITSHPTSSAQKIDVL KPRPRSFSDK NSCNDVNILP KVDVVPLHLG LNNNYVCESS ITCNKDEQKD KLININLLDG DNMWWESLLE ADVLGPEATE TAKGVTLPLD FEQIWARFDE ETELELN
Source	E.coli
Target Names	MYB113
Protein Names	Recommended name: Transcription factor MYB113 Alternative name(s): Myb-related protein 113 Short name= AtMYB113
Expression Region	1-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.