



Recombinant Arabidopsis thaliana Squamosa promoter-binding-like protein 11 (SPL11)

Product Code	CSB-EP861910DOA-B
Abbreviation	SPL11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FZK0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDCNMVSSSQ WDWEHLIMSN PSRTEDDSKQ LPTEWEIEKG EGIESIVPHF SGLERVSSGS ATSFWHTAVS KSSQSTSINS SSPEAKRCKL AESSPGDSC SNIDFVQVKA PTALEVSVAS AESDLCLKLG KRTYSEEYWG RNNNEISAVS MKLLTPSVVA GKSKLCGQSM PVPRCQIDGC ELDLSSAKGY HRKHKVCEKH SKCPKVS VSG LERRFCQQCS RFHAVSEFDE KKRSCRKRLS HHNARRRKPQ GVFSMNPERV YDRRQHTNML WNGVSLNARS EEMYEWGNNT YDTKPRQTEK SFTLSFQRGN GSEDQLVASS SRMFSTSQTS GGFPAGKSKF QLHGEDVGEY SGVLHESQDI HRALSLLSTS SDPLAQPHVQ PFSLLCSYDV VPK
Source	E.coli
Target Names	SPL11
Protein Names	Recommended name: Squamosa promoter-binding-like protein 11
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.